

MOLARITY

1. You are provided with the following

BA1 is a solution of sodium carbonate

BA2 is 0.2M hydrochloric acid solution.

In this experiment you are required to determine the concentration of sodium carbonate in moles per liter (molarity).

Procedure: pipette 20 or 25cm³ of sodium carbonate solution into a conical flask. Add 2drops of methyl orange indicator then titrate the mixture with BA2 from the Burette until the end point is reached.

Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used..... cm³

Final Burette reading cm ³			
Initial burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2 ,

Average volume of BA2

a) Find the number of moles of hydrochloric acid in BA2 that reacted.

.....

b) Write the equation for the reaction

.....

.....

c) Find the concentration of sodium carbonate in moles per liter (molarity)

.....

2. You are provided with the following

BA1 is sodium hydroxide solution

BA2 is 0.05M sulphuric acid solution.

In this experiment you are required to determine the concentration in moles per litre of sulphuric acid (molarity).

Procedure: pipette 20 or 25cm³ of BA1 in to a clean conical flask. Add 2drops of phenolphthalein indicator titrate the mixture with BA2 from the Burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results;

Volume of the pipette used..... cm³

Final burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

3. You are provided with the following

BA1 is 0.2M Sodium hydroxide solution.

BA2 is oxalic acid solution.

In this experiment you are required to determine the concentration of oxalic acid in the moles per liter (molarity)

Procedure; pipette 20 or 25cm³ of BA1 into a conical flask add 2 drops of phenolphthalein indicator titrate the mixture with BA2 from the burette until the end point is reached.

Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used..... cm³.

Final Burette reading cm ³			
Initial Burette readings cm ³			
Volume of BA2 used in cm ³			

Titre volumes of BA2,

Average volumes of BA2

a). Determine the number of moles of sodium hydroxide in BA1 that reacted

.....

.....

.....

.....

.....

.....

b).Write the equation for the reaction between BA1 and BA2

.....
.....
.....
.....
.....

c) .Determine the concentration of oxalic acid in moles per liter (molarity)

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....

4. You are provided with the following

BA1 is 0.1M sodium hydroxide solution.

BA2 is hydrochloric acid solution.

In this experiment you are required to determine the concentration in moles per litre of hydrochloric acid.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2

a) Determine the number of moles of BA1 that reacted.

.....

.....

.....

.....

b) Write the equation for the reaction between BA1 and BA2

.....

.....

.....

..... Determine the concentration of hydrochloric acid in BA2 in moles per litre (molarity)

.....

.....

.....

.....

.....

PERCENTAGE PURITY AND IMPURITY

1. You are provided with the following

BA1 is made by dissolving 12.7g of the impure sodium carbonate in one litre of solution.

BA2 is 0.2M Hydrochloric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium carbonate.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Find the number of moles of Hydrochloric acid in BA2 that reacted

.....

b) Write the equation for the reaction between BA1 and BA2

.....

c) Determine the concentration of sodium carbonate in moles per litre (molarity)

.....

.....

.....

.....

.....

..... Determine the percentage purity and impurity of sodium carbonate (Na=23, C=12, O=16)

.....

.....

.....

.....

.....

.....

2. You are provided with the following

BA1 is made by dissolving 6.4g of the impure sodium carbonate in 500cm³ of water to make a solution.

BA2 is 0.2M Hydrochloric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium carbonate.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used _____ cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Find the number of moles of Hydrochloric acid in BA2 that reacted

.....

.....
.....
.....
.....

b) Write the equation for the reaction between BA1 and BA2

.....
.....
.....

c) Determine the concentration of sodium carbonate in moles per litre (molarity)

.....
.....
.....
.....
.....
.....
.....
.....

d) Determine the percentage purity and impurity of sodium carbonate (Na=23, C=12, O=16)

.....
.....
.....
.....
.....
.....

3. You are provided with the following

BA1 is made by dissolving 3.175g of the impure sodium carbonate in 250cm³ of water to make a solution.

BA2 is 0.2M Hydrochloric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium carbonate.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used _____ cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....

Average volume of BA2.....

a) Find the number of moles of Hydrochloric acid in BA2 that reacted

.....

.....

.....

.....

.....

b) Write the equation for the reaction between BA1 and BA2

.....

.....

.....

c) Determine the concentration of sodium carbonate in moles per litre (molarity)

.....

.....

.....

.....

.....

.....

d) Determine the percentage purity and impurity of sodium carbonate (Na=23, C=12, O=16)

.....

4. You are provided with the following

BA1 is made by dissolving 10.2g of the impure Sodium hydroxide in one litre.

BA2 is 0.1M sulphuric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium hydroxide.

Procedure: Pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results

Volume of the pipette usedcm³

Final burette reading cm ³			
Initial burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2,

Average volume of BA2.....

a) Find the number of moles of BA2 that reacted

.....

.....
.....

b) Write the equation for the reaction

.....
.....
.....

c) Find the concentration of sodium hydroxide in moles per dm^3

.....
.....
.....
.....
.....
.....
.....
.....
.....

d) Find the percentage purity and percentage impurity of sodium hydroxide (Na=23, O=16, H=1)

.....
.....
.....
.....
.....
.....

5. You are provided with the following.

BA1 is made by dissolving 5.3g of the impure Sodium hydroxide in 500cm^3 of water to make a solution.

BA2 is 0.1M sulphuric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium hydroxide.

Procedure: Pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results

Volume of the pipette usedcm³

Final burette reading cm ³			
Initial burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2,

Average volume of BA2.....

e) Find the number of moles of BA2 that reacted

.....

.....

.....

.....

.....

f) Write the equation for the reaction

.....

.....

.....

g) Find the concentration of sodium hydroxide in moles per dm³

.....

.....

.....

.....

.....

.....

.....

.....

h) Find the percentage purity and percentage impurity of sodium hydroxide (Na=23, O=16, H=1)

.....

6. You are provided with the following

BA1 is made by dissolving 2.7g of the impure Sodium hydroxide in 250cm³ of water to make a solution.

BA2 is 0.1M sulphuric acid solution.

In this experiment, you are required to determine the percentage purity and impurity in sodium hydroxide.

Procedure: Pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results

Volume of the pipette usedcm³

Final burette reading cm ³			
Initial burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2,

Average volume of BA2.....

i) Find the number of moles of BA2 that reacted

.....

.....

j) Write the equation for the reaction

.....

.....

.....

k) Find the concentration of sodium hydroxide in moles per dm³

.....

.....

.....

.....

.....

.....

.....

.....

l) Find the percentage purity and percentage impurity of sodium hydroxide (Na=23, O=16, H=1)

.....

.....

.....

.....

.....

.....

7. You are provided with the following

BA1 is made by dissolving 13.4g of a mixture of sodium chloride and sodium carbonate in one litre.

BA2 is 0.2M hydrochloric solution.

In the experiment you are required to determine the percentage by mass of sodium chloride and percentage of sodium carbonate in the mixture.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of methyl orange indicator then

titrate the mixture with BA2 from the Burette until the end point is reached.

Record your results in the table below and repeat titration until you obtain the consistent results.

Volume of the reading used _____ cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2

Average volume of BA2

a) Find the number of moles of hydrochloric acid in BA2 that reacted.

.....

.....

.....

.....

.....

b) Write the equation for the reaction between BA1 and BA2.

.....

.....

.....

.....

c) Find the concentration of BA1 in moles per dm³ (molarity)

.....

.....

.....

.....

.....
.....
.....

d) Find the percentage by mass of sodium chloride and sodium carbonate in the mixture.

(Na=23, C=12, O=16)

.....
.....
.....
.....
.....
.....
.....
.....

8. You are provided with the following

BA1 is made by dissolving 13.6g of a mixture of Sodium sulphate and Sodium hydroxide in one litre.

BA2 is 0.05M sulphuric acid solution.

In this experiment, you are required to determine the percentage by mass of sodium sulphate and percentage of sodium hydroxide in the mixture.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of phenolphthalein indicator then titrate the mixture with BA2 from the Burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette usedcm³

Final Burette reading cm^3			
Initial Burette reading cm^3			
Volume of BA2 used cm^3			

Titre values of BA1,

Average volume of BA2

a) Find the number of moles of BA2 that reacted.

.....

.....

.....

.....

b) Write the equation for the reaction that took place.

.....

.....

.....

.....

c) Find the concentration of sodium – hydroxide in moles per liter (molarity)

.....

.....

.....

.....

.....

....._ d) Find the percentage by mass of sodium sulphate and percentage of sodium hydroxide in the mixture.

(Na=23, O=16, H=1).

.....

.....

.....

.....

.....

.....

.....

.....

.....

WATER OF CRYSTALLISATION

1. You are provided with the following

BA1 is made by dissolving 28.6g of hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ in one litre.

BA2 is 0.2M Hydrochloric acid. In this experiment, you are required to determine the value of n (water of crystallization) in hydrated sodium carbonate.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³				
Initial Burette reading cm ³				
Volume of BA2 used cm ³				

Titre values of BA2.....,

Average volume of BA2.....

d) Find the number of moles of BA2 that reacted

.....

.....

.....
.....
.....

e) Write the equation for the reaction between BA1 and BA2

.....
.....
.....

f) Determine the concentration of sodium carbonate in moles per litre

.....
.....
.....
.....
.....
.....
.....
.....
.....

g) Deduce the relative formula mass of Hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ and Hence find the value of n in the salt (Na=23, C=12, O=16, H=1)

.....
.....
.....
.....
.....

2. You are provided with the following

BA1 is made by dissolving 14.3g of hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ in 500cm^3 of water to make a solution.

BA2 is 0.2M Hydrochloric acid. In this experiment, you are required to determine the value of n (water of crystallization) in hydrated sodium carbonate.

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....

cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Find the number of moles of BA2 that reacted

.....

b) Write the equation for the reaction between BA1 and BA2

.....

c) Determine the concentration of sodium carbonate in moles per litre

.....

- d) Deduce the relative formula mass of Hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ and Hence find the value of n in the salt (Na=23, C=12, O=16, H=1)

.....

.....

.....

.....

.....

.....

3. You are provided with the following

BA1 is made by dissolving 7.15g of the hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ in 250cm^3 of water to make a solution.

BA2 is 0.2M Hydrochloric acid. In this experiment, you are required to determine the value of n (water of crystallization) in hydrated sodium carbonate.

Procedure: pipette 20 or 25cm^3 of BA1 into a conical flask, add 2 drops of methyl orange indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used..... cm^3

Final Burette reading cm^3			
Initial Burette reading cm^3			
Volume of BA2 used cm^3			

Titre values of BA2.....,

Average volume of BA2.....

- a) Find the number of moles of BA2 that reacted

.....

.....

.....
.....
.....

b) Write the equation for the reaction between BA1 and BA2

.....
.....

c) Determine the concentration of sodium carbonate in moles per litre

.....
.....
.....
.....
.....
.....
.....
.....

d) Deduce the relative formula mass of Hydrated sodium carbonate $\text{Na}_2\text{CO}_3 \cdot n\text{H}_2\text{O}$ and Hence find the value of n in the salt (Na=23, C=12, O=16, H=1)

.....
.....
.....
.....

4. You are provided with the following

BA1 is 0.2M sodium hydroxide solution.

BA2 is made by dissolving 12.6g of hydrated acid $\text{H}_2\text{CO}_4 \cdot X \text{H}_2\text{O}$ in one litre.

In this experiment you are required to determine the value of X (water of crystallization) in the acid

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Determine the number of moles of sodium hydroxide in BA1 used

.....

b) Write the equation for the reaction between BA1 and BA2

c) Determine the concentration of the acid in moles per dm³ (molarity)

d) Find the mass of one mole of the hydrated acid (molecular weight) and hence deduce the value of X in the acid. (H=1, C=12, O=16)

.....

.....

.....

.....

.....

.....

5 You are provided with the following

BA1 is 0.2M sodium hydroxide solution.

BA2 is made by dissolving 6.3g of hydrated acid $H_2C_2O_4 \cdot X H_2O$ in 500cm^3 of water to make a solution.

In this experiment you are required to determine the value of X (water of crystallization) in the acid

Procedure: pipette 20 or 25cm^3 of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used _____ cm^3

Final Burette reading cm^3			
Initial Burette reading cm^3			
Volume of BA2 used cm^3			

Titre values of BA2.....,

Average volume of BA2.....

e) Determine the number of moles of sodium hydroxide in BA1 used

.....

.....

.....

.....

.....
f) Write the equation for the reaction between BA1 and BA2

.....
.....
.....

g) Determine the concentration of the acid in moles per dm³ (molarity)

.....
.....
.....
.....
.....
.....
.....
.....
.....
.....

h) Find the mass of one mole of the hydrated acid (molecular weight) and hence deduce the value of X in the acid. (H=1,C=12,O=16)

.....
.....
.....
.....
.....
.....
.....
.....

6. You are provided with the following

BA1 is 0.2M sodium hydroxide solution. BA2 is made by dissolving 3.15g of hydrated acid $H_2C_2O_4 \cdot X H_2O$ in 250cm³ of water to make a solution..

In this experiment you are required to determine the value of X (water of crystallization) in the acid

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask, add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results

in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used _____ cm^3

Final Burette reading cm^3			
Initial Burette reading cm^3			
Volume of BA2 used cm^3			

Titre values of BA2.....,

Average volume of BA2.....

i) Determine the number of moles of sodium hydroxide in BA1 used

.....

j) Write the equation for the reaction between BA1 and BA2

.....

k) Determine the concentration of the acid in moles per dm^3 (molarity)

.....

l) Find the mass of one mole of the hydrated acid (molecular weight) and hence deduce the value of X in the acid. (H=1, C=12, O=16)

.....

MOLE RATIO

1. You are provided with the following

BA1 is 0.1M solution of base Q

BA2 is 0.05M solution of an acid P

In this experiment you are required to determine the mole ratio between Q and P. (stoichiometry).

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask and add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Find the number of moles of Q in BA1 that reacted

.....

.....

 b) find the number of moles of P in BA2

.....

 c) determine the mole ratio between Q and P

2. You are provided with the following

BA1 is 0.2M solution of base R

BA2 is made by dissolving 9.8g of an acid K in one litre.

In this experiment you are required to determine the mole ratio between Q and P. (stoichiometry).

Procedure: pipette 20 or 25cm³ of BA1 into a conical flask and add 2 drops of phenolphthalein indicator then titrate the mixture with BA2 from the burette until the end point is reached. Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette used.....cm³

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Titre values of BA2.....,

Average volume of BA2.....

a) Find the number of moles of R in BA1 that reacted

.....
.....
.....
.....
.....

b) find the number of moles of K in BA2

.....
.....
.....
.....
.....

c) determine the mole ratio between R and K

.....
.....
.....
.....
.....

BASICITY OF AN ACID

1. You are provided with the following
 BA1 is 0.1 M sodium hydroxide solution
 BA2 is 0.05M solution of an acid H_nX .

In this experiment you are required to determine the Basicity of an acid (value if n) in H_nX .

Procedure; Pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of phenolphthalein indicator then titrate it with BA2 from the Burette until the endpoint is reached, Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the Pipette

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Liter values of BA2,

Average volume of BA2

Find the number of moles of Sodium hydroxide in BA1 that reacted

.....

.....

b) Find the number of moles of the acid that reacted.

.....

c) Find the value of n in the acid.

.....

2. You are provided with the following

BA1 is 0.2 M sodium hydroxide solution

BA2 is 0.1M solution of an acid H_nX .

In this experiment you are required to determine the Basicity of an acid (value if n) in H_nX .

Procedure; pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of phenolphthalein indicator then titrate it with BA2 from the Burette until the endpoint is reached, Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the pipette.....

Final Burette reading cm ³			
Initial Burette reading cm ³			

Volume of BA2 used cm ³			
------------------------------------	--	--	--

Liter values of BA2,

Average volume of BA2

Find the number of moles of Sodium hydroxide in BA1 that reacted

.....

.....

.....

.....

.....

.....

b) Find the number of moles of the acid that reacted.

.....

.....

.....

.....

.....

c) Find the value of n in the acid.

.....

.....

.....

.....

3. You are provided with the following

BA1 is 0.3 M sodium hydroxide solution

BA2 is 0.2M solution of an acid H_nX .

In this experiment you are required to determine the Basicity of an acid (value if n) in H_nX .

Procedure; Pipette 20 or 25cm³ of BA1 into a conical flask add 2drops of phenolphthalein indicator then titrate it with BA2 from the Burette until the endpoint is reached, Record your results in the table below and repeat titration until you obtain consistent results.

Volume of the Pipette.....

Final Burette reading cm ³			
Initial Burette reading cm ³			
Volume of BA2 used cm ³			

Liter values of BA2,

Average volume of BA2

Find the number of moles of Sodium hydroxide in BAA1 that reacted

.....

b) Find the number of moles of the acid that reacted.

.....

c) Find the value of n in the acid.

.....

.....

.....

.....

.....

THERMO CHEMISTRY

Determination of molar enthalpy of neutralization of hydrochloric acid by NaOH solution

NOTE: molar enthalpy of neutralization is the heat change which occurs when one mole of water is formed. When one mole of hydroxide ion from an alkali reacts with one mole of hydrogen ion from an acid.

You are provided with solutions

BA1 which is 0.05M sodium hydroxide solution

BA2 which is 0.5M hydrochloric acid solution

Method

- Pipette 20.0 or 25.0cm³ of BA1 into a conical flask
- Add 5cm³ of hydrochloric acid from the burette, stir the mixture using a thermometer and note the final temperature of the mixture
- Repeat the steps 1 and 2 using 10, 15, 20, 25, 30 and 35cm³ of hydrochloric acid
- Record your results in the table below including the initial temperature of hydrochloric acid.

Results

Volume of pipette used25.0cm³.....

Initial temperature of HCl 23.0°C

Volume of HCl added(cm ³)	0	5	10	15	20	25	30	35
Volume of NaOH added (cm ³)	25	25	25	25	25	25	25	25
Final temperature of mixture(°C)	23.3	24.0	24.6	25.3	25.8	26.4	25.7	25.2

- i) From the table, determine the volume of HCl required to completely neutralize to NaOH, state the temperature at the point of complete neutralize 25.00cm³
- ii) Calculate the temperature rise at the point of neutralization
- iii) In the table, before adding any acid, i.e. when volume of HCl is zero, the temperature recorded is the initial temperature of NaOH
- iv) Calculate the heat gained by the mixture
- v) Determine the moles of HCL that reacted
- vi) Determine the molar enthalpy of neutralization of HCl

RATES OF REACTIONS

To investigate the effect of concentration on the rate of reaction between dilute hydrochloric acid and sodium thiosulphate solution

Requirements

2M HCl in a burette

0.15M sodium thiosulphate

Stop clock

Procedure

- i) Mark a large cross (x) on a piece of white paper
- ii) Pour 50cm³ of 2M HCl from the burette into the beakers containing sodium thiosulphate solution and immediately start a stop clock
- iii) Look through the solution in the flask at the cross and note the time taken for the cross to disappear
- iv) Repeat the experiment using different volumes of thiosulphate as shown and making sure that the volumes are adjusted with distilled water.

Table

Volume of sodium thiosulphate cm ³	Volume of water added	Volume of acid	Time t/s	1/t (s ⁻¹)
50	0	50		
40	10	50		
30	20	50		
20	30	50		
10	40	50		

Question

- a) Plot a graph of volume of thiosulphate (y-axis) against $\frac{1}{t}$ (x-axis)
- b) Use your graph to explain how the rate of reaction depends on the concentration of $S_2O_3^{2-}$
- c) In this experiment $\frac{1}{t}$ is used as a measure of the rate of reaction. Explain how.

NB: 0.15M $Na_2S_2O_3$ is prepared by dissolving 37.2g of the salt in distilled water and the solution made to 1000cm³

QUALITATIVE ANALYSIS

Qualitative analysis is the systematic identification of chemical constituents of a particular substance in form of ions (cations and anions)


Cations carry positive charges

Anions carry negative charges


COMMON CATIONS IDENTIFIED

- a) Transition metal ions


These appear as coloured Cations and include:-



- b) Non-coloured Cations



COMMON ANIONS IDENTIFIED



PRELIMINARY TEST

SOLIDS

Table 1

OBSERVATION	INFERENCE
Black colour	An oxide or sulphide
Yellow colour	Lead(II) oxide or iron(III) salt
Green colour	Iron(II) or copper(II) salt

Blue colour	Iron (//) or copper (//) salt
Smell of ammonia	Ammonium salt
Smell of hydrogen sulphide	Sulphide
Smell of sulphur dioxide	sulphite
deliquescent	Chloride or nitrate

FLAME TESTS

Place little of the substance on a watch glass, moisten it with pure concentrated hydrochloric acid, and heat little on a clean platinum or alchrome wire. Note the colour of the flame if you think potassium is present at flame through glass.

Table 2

OBSERVATION	INFERENCE
Brilliant yellow	Sodium
Red	Calcium
Blue	Lead
Blue-green	Copper
Lilac (crimson thru blue glass)	Potassium

ACTION OF HEAT ON SOLID

Table 3

	OBSERVATION	DEDUCTION
I)	Colourless liquid which turns anhydrous copper (II) sulphate blue forms on cooler parts of the tube	Solid contains water of crystallization of hydrated salt
II)	White sublimate forms	NH_4^+ present
III)	Cracking sound is made	Pb^{2+}
IV)	Residue yellow when hot and white on cooling	ZnO
V)	Residue yellow when hot and colour persists on cooling	PbO

VI)	Residue black	CuO, FeO, NiO suspected
-----	---------------	-------------------------

GAS**Table 4**

GAS	COLOUR	SMELL	ACTION WITH SPLINT	LITMUS
Chlorine	Greenish Yellow	Irritating characteristic	Doesn't burn or support combustion	Turn Red then is bleached
Hydrogen chloride or (bromine)	Colourless misty fumes in damp air	Irritating	Doesn't burn or support combustion	Red
Nitrogen dioxide	Reddish brown	Irritating characteristic	Doesn't burn or support combustion	Red
Nitric acid vapour	Pale yellow fumes	Irritating	Split burns	Red
Sulphur dioxide	White fumes	Irritating	Doesn't burn or support	Red

GAS	COLOUR	SMELL	ACTION WITH SPLINT	LITMUS
Hydrogen sulphide	None	Rotten eggs	Burns with blue flame, depositing yellow sulphur on cold surface	Faint red
Ammonia	None	Chocking characteristics	Doesn't burn or support combustion	Blue
Hydrogen	None	None(if pure)	Burns with blue	None

			flame if mixed with air	
Oxygen	None	None	Relights glowing splint	None
Carbon dioxide	None	Faint (not easy to detect)	Doesn't burn or support combustion	Red
Water vapour	None	None	Doesn't burn or support combustion	None
Nitrogen	None	None	Doesn't burn or support combustion	None

NB: You are advised to test for gas whenever;

- i) You heat the solid
- ii) You add sodium hydroxide solution to sample and heat
- iii) You add acid to the sample

DISSOLVING A SOLID IN WATER

It is important to observe whether the substance is soluble, partially soluble or insoluble in water.

If the substance is insoluble or partially soluble then filtering is done.

Soluble	Insoluble or partially soluble
All nitrates	
All common chlorides	AgCl and PbCl ₂ (PbCl ₂ is insoluble in hot water or on warming)
All common sulphates	BaSO ₄ , PbSO ₄ , CaSO ₄ is sparingly soluble
K ₂ CO ₃ , Na ₂ CO ₃ , (NH ₄) ₂ CO ₃	All common carbonates are insoluble.

Example

Observation	Deduction
Solid dissolves forming a colourless solution	Al ³⁺ , Zn ²⁺ , Mg ²⁺ , NH ₄ ⁺
Solid dissolves forming a green or pale	Fe ²⁺ , Cu ²⁺ present

green solution	
Solid dissolves forming a blue solution	Cu^{2+} present

DISSOLVING A SOLID IN WATER AND FILTER

In this case, you state the colour of residue and filtrate

a) Filtrate

Table 5

OBSERVATION	DEDUCTION
Filtrate is colourless (not colourless solution or suspension)	Al^{3+} , Zn^{2+} , Pb^{2+} , Mg^{2+} , Ca^{2+} or Ba^{2+} present Ba^{2+} , present
Filtrate is blue	Cu^{2+} present
Filtrate is green	Cu^{2+} , Fe^{2+} present
Filtrate green, turn brown/yellow at surface on standing	Fe^{2+} , present, oxidized to Fe^{3+}
Filtrate is brown/yellow	Fe^{3+} , present

Residue

Table 6

OBSERVATION	DEDUCTION
Residue is white (not white solution or suspension)	Al^{3+} , Zn^{2+} , Pb^{2+} , Mg^{2+} present
Residue is blue	Cu^{2+} present
Residue is green	Cu^{2+} , Fe^{2+} present
Residue green, turn brown/yellow at surface on standing	Fe^{2+} , present, oxidized to Fe^{3+}
Residue is brown/yellow	Fe^{3+} , present

REACTION WITH SODIUM CARBONATE SOLUTION**Table 7**

OBSERVATION	DEDUCTION
No observable change	Pb^{2+} , present
White precipitate forms and a colourless acidic gas which turns lime water milky is given off	The gas is CO_2 Therefore Al^{3+} confirmed present
A colourless gas which turns red litmus paper blue is given off (on warming or cold) NB. Here no ppt is formed	The gas is NH_3 Therefore NH_4^+ confirmed present
White precipitate is formed	Zn^{2+} , Mg^{2+} , (Ca^{2+} , Ba^{2+}) present
Blue/ dirty blue precipitate turns black on heating	Cu^{2+} confirmed present
Red brown/ rust brown precipitate forms and a colourless acidic gas which turns lime water milky is given off	The gas is CO_2 Therefore Fe^{2+} , confirmed present
Green precipitate forms, insoluble in excess (No effervescence)	Fe^{2+} confirmed present

ACTION OF SODIUM HYDROXIDE SOLUTION


This tests for cations only

Table 8

OBSERVATION	DEDUCTION
White precipitate, insoluble in excess	Ca^{2+} , Mg^{2+} present
White precipitate soluble in excess forming colourless solution	Pb^{2+} , Zn^{2+} , Al^{3+} present
Blue precipitate, insoluble in excess, or ppt turns black on heating	Cu^{2+} present
Green precipitate, insoluble in excess, or	Fe^{2+} , present, oxidized to Fe^{3+}

ppt turns brown on adding hydrogen peroxide, turns brown on standing	
Dark brown precipitate, insoluble in excess	Fe^{3+} , present
No observable change, on warming a colourless gas which turns red litmus blue evolved (ammonia gas)	NH_4^+ , present

NB: Hydroxides of lead, zinc and aluminium dissolve in excess sodium hydroxide solution.



ACTIONS WITH AMMONIA SOLUTION

Note: Ammonia solution test for only cations

Table 9

OBSERVATION	DEDUCTION
White precipitate formed, insoluble in excess	Al^{3+} , Pb^{2+} , Mg^{2+} Confirmed present
White precipitate forms soluble in excess	Zn^{2+} , Confirmed present

Pale blue / blue precipitate formed soluble in excess alkali, forming a deep blue solution	Cu^{2+} , Confirmed present
Dirty/dark green precipitate formed insoluble in excess alkali, precipitate turns brown on standing	Fe^{2+} , confirmed present being oxidized to Fe^{3+}
Rust brown precipitate forms, insoluble in excess ammonia solution	Fe^{3+} , Confirmed present

Copper (//) Hydroxide dissolve in excess aqueous ammonia to complex compounds.





Tetra-amine copper (//) and tetra-amine zinc ions

ACTION OF DILUTE HYDROCHLORIC ACID

Add cold dilute hydrochloric acid to the solid substance in a test tube. If there is no reaction, warm gently.

Table 10

OBSERVATION	INFERENCE
Colourless acidic gas which turns lime water milky evolved	The gas is CO_2 CO_3^{2-} confirmed present
Sulphur dioxide	sulphite
Hydrogen sulphide	Sulphide



ADDITION OF LEAD NITRATE OR ETHANOATE SOLUTION

Note: These reagents confirm the presence of SO_4^{2-} , CO_3^{2-} / $\text{C}_2\text{O}_4^{2-}$ ions

Table 11

OBSERVATION	DEDUCTION
White precipitate forms, dissolves on heating and re-precipitate on cooling	Cl^- ion confirmed present
White precipitate forms (no effect on heating), which dissolves in nitric acid forming a colourless solution with effervescence of a colourless acidic gas which turns lime water milky	The gas is CO_2 CO_3^{2-} ions
White precipitate forms (no effect on heating) which is insoluble in nitric acid	SO_4^{2-} confirmed present

REACTION WITH POTASSIUM IODIDE SOLUTION**Table 12**

OBSERVATION	DEDUCTION
Yellow precipitate formed	Pb^{2+} confirmed present
Brown precipitate forms (form blue solution)	Cu^{2+} confirmed present

CONFIRMATORY TESTS FOR CATIONS

TEST	OBSERVATION	CATION CONFIRMED PRESENT
Add 1cm ³ of dilute NaOH solution and warm	A colourless gas which turns red litmus paper blue is given off	NH_4^+ present
Add NH_3 solution drop wise until in excess	White precipitate forms which dissolves in excess alkali forming a colourless solution	Zn^{2+} present
Add few drops of Potassium hexacyanoferrate (II) solution	White precipitate forms	Zn^{2+} present

Add 2-3 drops of Na_2CO_3 solution	White precipitate forms	Zn^{2+} present
Add 1cm ³ of dilute HCl acid	White precipitate soluble on heating and re-appears on cooling	Pb^{2+} present
Add 3 drops of Potassium Iodide Solution	Yellow precipitate forms	Pb^{2+} present
Add dilute Sodium hydrogen phosphate solution	White gelatinous precipitate forms	Al^{3+} present
Add dilute ammonium solution and ammonium chloride solution followed by a little dilute Sodium Hydrogen phosphate solution	A white crystalline precipitate forms	Mg^{2+} present

SUMMARY OF CONFIRMATORY TESTS FOR Fe^{2+} , Fe^{3+} AND Cu^{2+}

	Fe^{2+}	Fe^{3+}	Cu^{2+}
Add 3 drops of Potassium hexacyano ferrate(II) solution	Dark blue precipitate	dark blue precipitate	dark brown precipitate
Add Potassium/ sodium/ ammonium thiocyanate solution		Deep red/ blood red solution forms	
Add sodium carbonate solution	Dirty green precipitate insoluble in excess alkali	Rust brown precipitate insoluble in excess in alkali	Blue precipitate soluble in excess alkali forming deep blue solution
add sodium carbonate	green precipitate	brown precipitate	Blue precipitate

solution		with effervescence of CO ₂ gas	turns black on heating
----------	--	---	------------------------

CONFIRMATORY ANALYSIS TEST FOR ANIONS

ION	TEST	OBSERVATION
SO ₃ ²⁻	To the solution add dilute HCl and warm	Colourless gas evolved which turns acidified potassium dichromate solution from yellow to green
CO ₃ ²⁻	Add HCl or HNO ₃ or H ₂ SO ₄	Effervescence of a colourless gas that turns lime water milky
NO ₃ ⁻	Add concentrated Sulphuric acid and cold freshly prepared aqueous Iron (II) sulphate solution To the cold solution add a few Copper turnings and then concentrated Sulphuric acid and warm	Brown ring formed at the interface Blue solution forms and brown fumes are given off
Cl ⁻	To the solution add dilute Nitric acid and then Silver Nitrate solution followed by dilute Ammonia solution and shake	A white ppt is formed insoluble in the acid. Which dissolves aqueous ammonia forming a colourless solution
SO ₄ ²⁻	To the solution add few drops acidified aqueous Barium Nitrate. To the	White ppt is formed insoluble in the acid

	solution add dilute Hydrochloric acid then Barium chloride solution	White ppt insoluble in the acid
--	---	---------------------------------

HOW TO DIFFERENTIATE BETWEEN CO_3^{2-} AND HCO_3^-

To the solution add few solids of Magnesium Sulphate	<p>a) White ppt soluble indicates that HCO_3^- ions are present</p> <p>b) White ppt insoluble indicates that CO_3^{2-} ions are present</p>
--	---

Table 13

TEST	OBSERVATION	DEDUCTION
Add dilute HCl acid, HNO_3 or H_2SO_4	Effervescence of a colourless acidic gas which turns lime water milky	The gas is CO_2 CO_3^{2-} confirmed present
Add barium chloride then dilute HNO_3	White precipitate soluble in HNO_3 with Effervescence of a colourless acidic gas that turns lime water milky	The gas is CO_2 CO_3^{2-} confirmed present
Add barium chloride then dilute HCl acid	White precipitate soluble in HNO_3 with Effervescence of a colourless acidic gas that turns lime water milky	The gas is CO_2 CO_3^{2-} confirmed present
Add lead nitrate or lead	White precipitate soluble in	The gas is CO_2

ethanoate solution, then dilute HNO_3	HNO_3 with Effervescence of a colourless acidic gas that turns lime water milky	CO_3^{2-} confirmed present
Add lead nitrate or or lead ethanoate solution, then dilute nitric acid	White precipitate insoluble in nitric acid	SO_4^{2-} confirmed present
Add lead nitrate solution then dilute hydrochloric acid	White precipitate insoluble in hydrochloric acid	SO_4^{2-} confirmed present
Add dilute nitric acid then a few drops of silver nitrate solution followed by dilute ammonium solution hydrochloric acid	White precipitate insoluble in nitric acid	Cl^- confirmed present
To a solid sample, add concentrated H_2SO_4 and heat	Effervescence occurs with misty fumes which turn blue litmus paper red and form dense white fumes which are ammonia	Cl^- confirmed present

TEST	OBSERVATION	DEDUCTION
Add equal volume of freshly prepared iron(II) sulphate solution then the tube is tilted and concentrated sulphuric acid added down the side of the tubes	Two separate layers form a brown ring appear between them	NO_3^- confirmed present
Heat the solid sample with concentrated sulphuric acid	Brown fumes (of nitrogen dioxide) forms	NO_3^- confirmed present
Add magnesium sulphate solution (or magnesium chloride solution and boil)	No observable change but on boiling a white precipitate forms. NB: CO_3^{2-} forms white in cold	HCO_3^- confirmed present

1. You are provided with substance Q which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas (es) evolved. Record your observations and deduction in the table below.

Test	Observation	Deduction
a) Heat a spatula endful of Q in a dry test tube to constant mass. Or until there is no further change.		
b) To another spatula endful of Q put in to a test tube add about 5cm ³ of distilled water, shake well, filter, and keep both the residue and the filtration.		

<p>c) Divide the filtrate into five parts</p> <p>i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		
<p>iii) To the third part add lead (II) nitrate solution</p>		
<p>iv) To the fourth part add lead(II) nitrate solution followed by dilute nitric acid</p>		

<p>v) To the fifth part carry out a test of your own choice to identify the Anions in the filtration.</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		
<p>d) To the residue put in to the test tube add dilute nitric acid until there is no further change</p> <p>Then divide the resultant solution in to two parts</p> <p>i) To the first part add sodium hydroxide solution drop wise until in excess.</p>		
<p>d)ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		

--	--	--

Cations in Q,

Anions in Q,

2) You are provided with substances X which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

Test	Observation	Deduction
------	-------------	-----------

<p>a) Heat a spatula end full of X in a dry test tube to constant mass or until there is no further change.</p>		
<p>b) To another spatula endful of X put in to a test tube add about 5cm³ of distilled water, shake well, filter, and keep both the residue and the filtration.</p>		
<p>c) Divide the filtrate into five parts i) To the first part add sodium hydroxide solution drop wise until in excess</p>		

<p>ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		
<p>iii) To the third part add lead(II) nitrate solution</p>		
<p>(iv) To the fourth part add lead (II) nitrate solution followed by dilute nitric acid.</p>		
<p>(v) To the fifth part carry out a test of your choice to identify the Anions in the filtration</p> <p>.....</p> <p>.....</p> <p>.....</p>		

<p>.....</p> <p>.....</p> <p>.....</p>		
<p>(d) To the residue put in to the test tube add dilute nitric acid until there is no further change. Then divide the resultant solution in two parts.</p> <p>i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>d) To the second part add aqueous ammonia solution drop wise until in excess</p>		

Cations in X,

Anions in X

3) You are provided with substances Y which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

Test	Observation	Deduction
a) Heat a spatula end full of Y in a dry test tube to constant mass or until there is no further change.		
b) To another spatula endful of Y put in to a test tube add about 5cm ³ of distilled water, shake well, filter, and keep both the residue and the filtration.		

<p>c) Divide the filtrate into five parts i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		
<p>iii) To the third part add lead(II) nitrate solution</p>		

<p>(iv) To the fourth part add lead (II) nitrate solution followed by dilute nitric acid.</p>		
<p>(v) To the fifth part carry out a test of your choice to identify the Anions in the filtration</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		
<p>(d) To the residue put in to the test tube add dilute nitric acid until there is no further change. Then divide the resultant solution in two parts.</p> <p>i) To the first part add sodium hydroxide solution drop wise until in excess</p>		

d)(II) To the second part add aqueous ammonia solution drop wise until in excess		
--	--	--

Cations in Y

Anions in Y

4) You are provided with substances W which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

Test	Observation	Deduction
------	-------------	-----------

<p>a) Heat a spatula end full of W in a dry test tube to constant mass or until there is no further change.</p>		
<p>b) To another spatula endful of W put in to a test tube add about 5cm³ of distilled water, shake well, filter, and keep both the residue and the filtration.</p>		
<p>c) Divide the filtrate into five parts i) To the first part add sodium hydroxide solution drop wise until in excess</p>		

<p>ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		
<p>iii) To the third part add lead(II) nitrate solution</p>		
<p>(iv) To the fourth part add lead (II) nitrate solution followed by dilute nitric acid.</p>		
<p>(v) To the fifth part carry out a test of your choice to identify the Anions in the filtration</p> <p>.....</p> <p>.....</p> <p>.....</p>		

.....		
(d) To the residue put in to the test tube add dilute nitric acid until there is no further change. Then divide the resultant solution in two parts. i) To the first part add sodium hydroxide solution drop wise until in excess		
d)(II) To the second part add aqueous ammonia solution drop wise until in excess		

Cations in W,

Anions in W.....

5) You are provided with substances P which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

Test	Observation	Deduction
------	-------------	-----------

<p>a) Heat a spatula end full of P in a dry test tube to constant mass or until there is no further change.</p>		
<p>b) To another spatula endful of P put in to a test tube add about 5cm³of distilled water, shake well, filter, and keep both the residue and the filtration.</p>		
<p>c) Divide the filtrate into five parts i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess.</p>		

<p>iii) To the third part add lead(II) nitrate solution</p>		
<p>(iv) To the fourth part add lead (II) nitrate solution followed by dilute nitric acid.</p>		
<p>(v) To the fifth part carry out a test of your choice to identify the Anions in the filtration</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		
<p>(d) To the residue put in to the test tube add dilute nitric acid until there is no further change. Then divide the resultant solution in two parts.</p> <p>i) To the first part add sodium hydroxide solution dropwise until in excess</p>		

d)(II) To the second part add aqueous ammonia solution drop wise until in excess		
--	--	--

Cations in P,

Anions in P

6) You are provided with substance A which contains two cations and two anions. You are required to carry out the given tests to identify them any gase(s) evolved. Record your observation and deductions in the table below.

Test	Observation	Deduction
a) Heat a spatula endful of A in a dry test tube to constant mass		
b) To another spatula endful of A put into a test tube add about 5cm ³ of		

di distilled water shake well, filter, and keep both the residue and filtrate.		
c) Divide the filtration into four parts i) To the first part add sodium hydroxide solution drop wise until in excess.		
ii) To the second part add ammonia solution drop wise until in excess.		
iii) To the third part add lead (II) nitrate solution followed by dilute nitric acid.		
iv) To the fourth part carryout your own test of your choice to identify the anion in the filtration.		

<p>d) To the residue put into the test tube add dilute nitric acid solution until there is no further change. Divide the resultant solution into three parts.</p> <p>(i) To the first part add sodium hydroxide solution dropwise until in excess.</p>		
<p>(ii) To the second part add aqueous ammonia solution dropwise until in excess.</p>		
<p>(iii) To the third part carry out your own test to identify the cations in the residue.</p>		

Cations in A,

Anions in A

7) You are provided with substance K which contains TWO cations and two Anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

Test	Observation	Deduction
------	-------------	-----------

<p>a) Heat spatula endful of K in a dry test tube to constant mass.</p>		
<p>b) To another spatula endful of K put into a test tube add about 5cm³ of water shake well, filter keep both the residue and the filtration</p>		
<p>c) Divide the filtrate into four parts add sodium hydroxide solution and heat.</p>		

<p>ii) To the second part add lead(II) nitrate solution followed by dilute nitric acid solution</p>		
<p>iii) To the third part add lead (II) nitrate solution and warm the mixture.</p>		
<p>v) To the fourth part carryout your own test to identify the anion in the filtration.</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		

.....		
<p>d) To the residue put into a test tube add dilute nitric acid solution and divide the resultant solution into three parts</p> <p>i) To the first part add sodium hydroxide solution dropwise until in excess</p>		
<p>ii) To the second part add aqueous ammonia solution dropwise until in excess</p>		

<p>iii) To the third part carry out your own test to identify the cation in the residue.</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		
---	--	--

Cations in K.....,

Anions in K

8) You are provided with substance R which contains two cations and two anions. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observation and deductions in the table below.

Test	Observation	Deduction
<p>a) heat a spatula endful of R in a dry test tube to constant mass</p>		

<p>b) To another spatula endful of R put into a test tube add about 5cm³ of water, shake well, Filter and keep both the residue and the filtration</p>		
<p>c) Divide the filtration into four parts</p> <p>i) add sodium hydroxide solution and heat</p>		
<p>ii) to the second part add lead nitrate solution followed by dilute nitric acid</p>		
<p>iii) To the third part add lead nitrate solution and heat</p>		
<p>iv) To the fourth part carry out your own test to identify the anions in the filtration</p>		

.....		
Put the residue in the test tube add dilute nitric acid solution then Divide the resultant solution into two parts To the first part add sodium hydroxide solution drop wise until in excess		
ii) To the second part add aqueous ammonia solution drop wise until in excess		

Cations in R..... ,

Anions in R.....

9. You are provided with substance Q which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observation	Reduction
------	-------------	-----------

a) Heat a spatula endful of Q in a dry test tube to constant mass		
b) To another spatula endful of Q put into a test tube add about 5cm ³ of water. Shake well to make a solution		
c) Divide the resultant solution into five parts i) Add sodium hydroxide solution drop wise until in excess and warm		
ii) To the second part add aqueous ammonia solution drop wise until in excess		
iii) To the third part carry out your own test to identify one of the cation		

iv) To the fourth part add lead (//) nitrate solution followed by dilute nitric acid		
v) To the fifth part carry out your own test to identify the anion in the filtrate		

10. You are provided with substance T which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observat i on	Reduct i on
a) Heat a spatula end ful of T in a dry test tube to constant mass		

<p>b) To another spatula endful of T put into a test tube add about 5cm³ of water. Shake well to make a solution</p>		
<p>c) Divide the resultant solution into five parts</p> <p>i) Add sodium hydroxide solution drop wise until in excess and heat</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess</p>		
<p>iii) To the third part add lead nitrate solution followed by dilute nitric acid</p>		

iv) To the fourth part carry out your own test to identify one of the cation		
v) To the fifth part carry out your own test to identify the anion in the filtrate 		

Cations in T..... ,

Anions in T

11. You are provided with substance N which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observati on	Reducti on
------	--------------	------------

<p>a) Heat a spatula endful of N in a dry test tube to constant mass</p>		
<p>b) To another spatula endful of N put into a test tube add about 5cm³ of water. Shake well to make a solution</p>		
<p>c) Divide the resultant solution into five parts</p> <p>i) Add sodium hydroxide solution drop wise until in excess and heat</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess</p>		

iii) To the third part add lead nitrate solution followed by dilute nitric acid		
iv) To the fourth part carry out your own test to identify one of the cation		
v) To the fifth part carry out your own test to identify the anion in the filtrate		

Cations in N..... ,

Anions in N,

12. You are provided with substance E which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observation	Reduction
------	-------------	-----------

<p>d) Heat a spatula endful of E in a dry test tube to constant mass</p>		
<p>e) To another spatula endful of E put into a test tube add about 5cm³ of water. Shake well to make a solution</p>		
<p>f) Divide the resultant solution into five parts</p> <p>vi) Add sodium hydroxide solution drop wise until in excess and heat</p>		
<p>vii) To the second part add aqueous ammonia solution drop wise until in excess</p>		

viii) To the third part add lead nitrate solution followed by dilute nitric acid		
ix) To the fourth part carry out your own test to identify one of the cation		
x) To the fifth part carry out your own test to identify the anion in the filtrate		

Cations in E..... ,

Anions in E

13. You are provided with substance D which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observat i on	Reduct i on
------	---------------	-------------

<p>a) Heat a spatula endful of D in a dry test tube to constant mass</p>		
<p>b) To another spatula endful of D put into a test tube add about 5cm³ of water. Shake well to make a solution</p>		
<p>c) Divide the resultant solution into five parts</p> <p>i) Add sodium hydroxide solution drop wise until in excess and heat</p>		
<p>To the second part add aqueous ammonia solution drop wise until in excess</p>		

ii) To the third part add lead nitrate solution followed by dilute nitric acid		
iii) To the fourth part carry out your own test to identify one of the cation		
iv) To the fifth part carry out your own test to identify the anion in the filtrate		

Cations in D..... ,

Anions in D.....,

14. You are provided with substance X which contains two Cations and one anion. You are required to carry out the given tests to identify them and any gase(s) evolved, Record your observations and deductions in the table below

Test	observat i on	Reduct i on
a) Heat a spatula end ful of X in a dry test tube to constant mass		

<p>b) To another spatula endful of X put into a test tube add dilute nitric acid and warm until there is no further change.</p> <p>Then add sodium hydroxide solution dropwise until in excess.</p> <p>Filter, keep both the residue and the filtrate</p>		
<p>c) To the filtrate add dilute nitric acid until solution is just acidic then divide it into three parts</p> <p>i) To the first part of the acidified filtrate add sodium hydroxide solution dropwise until in excess</p>		
<p>ii) To the second part of the acidified filtrate add ammonia solution dropwise until in excess</p>		
<p>iii) To the third part of the acidified filtrate add, carry out your own test to identify the cations in the filtrate</p> <p>.....</p> <p>.....</p> <p>.....</p>		

.....		
<p>d) To the residue put into a test tube then dissolve it into dilute nitric acid solution</p> <p>Then divide the resultant solution into three parts.</p> <p>i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>ii) To the second part add iron filling shake well and allow the mixture to stand for about 3 minutes</p>		
<p>iii) To the third part add aqueous ammonia solution drop wise until in excess.</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		

Cations in X..... ,

Anions in X.....

15) You are provided with substances Y which contain two cations and one anions. Carry out the given tests to identify the cations and anions in Y. Record your observations and deductions in the table below.

Test	Observation	Deduction
a) Heat a spatula end full of Y in a dry test tube to constant mass		
b) To the residue in (a) above allow to cool then add dilute nitric acid and warm until there is no further exchange then add sodium hydroxide solution drop wise until in excess. Filter , keep both the residue and filtrate		
c) To the filtration add dilute nitric acid until the solution is just acidic then divide it into four parts. i) To the first part add sodium hydroxide solution drop wise until in excess.		

<p>ci) To the second part of acidified filtrate add ammonia solution drop wise until in excess</p> <p>iii To the third part add lead nitrate solution</p>		
<p>iv) To the fourth part add Barium nitrate solution</p>		
<p>d) Dissolve the residue into dilute nitric acid solution. Then divide it into three parts</p> <p>i) To the first part add zinc powder and shake then allow the mixture to stand for about 2 minutes</p>		
<p>ii) To the second part add sodium hydroxide solution drop wise until in excess</p>		

iv) To the third part add ammonia solution drop wise in excess.		

Cations in y..... ,

Anions in y.....

16. You are provided with substance F which contains two cations and one anion. You are required to carry out the given tests to identify them and any gas(es) evolved. Record your observations and deductions in the table below.

a) Heat a spatula endful of F in a dry test tube until there is no further change		
b) To the another spatula endful of F		

<p>put into a test tube add dilute nitric acid solution until there is no further change</p> <p>then add sodium hydroxide solution drop wise until in excess filter, keep both the residue and the filtration</p>		
<p>c) To the filtration add dilute nitric acid in the solution then divide it into three parts</p> <p>i) to the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>iii) To the third part carry out your own test to identify the cations in the filtration.</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p> <p>.....</p>		

<p>d) To the residue put it into a test tube add dilute nitric acid to make a solution Divide the resultant solution into two parts</p>		
<p>i) To the first part add sodium hydroxide solution drop wise until in excess</p>		
<p>ii) To the second part add aqueous ammonia solution drop wise until in excess</p>		

Cations in F..... ,

Anions in F.....